

TORNEIRAS VALVULADAS

AUTOLIMPANTES



GB-04008 Torneira Bi-Valvulada Alto Fluxo



GB-06005B Extensor 6Fr x 5cm c/ Torneira Bi-Valvulada



GB-05018 Extensor 12Fr x 10cm c/ Torneira Bi-Valvulada



GB-08120V Extensor 8Fr x 120cm c/ Torneira Valvulada

GB-12060V Extensor 12Fr x 60cm c/ Torneira Valvulada

GB-12120V Extensor 12Fr x 120cm c/ Torneira Valvulada



GB-12020B Extensor 12Fr x 20cm c/ Torn. Bi-Valv. Alto Fluxo

GB-12060B Extensor 12Fr x 60cm c/ Torn. Bi-Valv. Alto Fluxo

GB-12120B Extensor 12Fr x 120cm c/ Torn. Bi-Valv. Alto Fluxo



GB-04009 Torneira Valvulada



GB-04028 Extensor 12Fr x 20cm c/ 2 Torneiras Valvuladas



GB-04038 Extensor 12Fr x 20cm c/ 3 Torneiras Valvuladas



GB-04138 Extensor 8Fr x 120cm c/ 3 Torneiras Valvuladas



GB-04158 Extensor 8Fr x 120cm c/ 5 Torneiras Valvuladas


TORNEIRAS VALVULADAS AUTOLIMPANTES

Alto Fluxo Gravitacional
aproximado de 12L/h


NOVO
rotor de
segurança



Exclusivo sistema fechado de passagem de fluidos que permite limpeza interna automática quando em fluxo contínuo, garantindo mínimo volume residual.



Volume Residual



Volume Residual

